HIRSCHMANN®

Technical data sheet

Article: Graduated pipette 1 ml, Class AS
Code no.: 110 01 33

Material:
AR (soda lime glass)

Chemical properties:

- Hydrolytic resistance at 98°C: Class 3 acc. ISO 719
- Hydrolytic resistance at 121°C: Class 3 acc. ISO 720
- Acid resistance: acc. ISO 1776
- Alkali resistance: Class 2 acc. ISO 695

Class:
AS

Hirschmann Laborgeräte
GmbH & Co.KG

Roland Benz
Quality Manager

Eberstadt, 25.07.2022