Dimethyl Sulfoxide
Baker Analyzed

Material No: 7033.5000JC
Effective from: 2013/07/22

Product Specification

Test	Specification
Assay	>= 99.0 %
Color (APHA)	<= 10
Density (g/ml) at 25°C	>= 1.095
Freezing Point	18.0-19.0 °C
Residue after Evaporation	<= 0.01 %
Water (H2O)	<= 0.2 %

For Laboratory, Research or Manufacturing use.

Storage conditions: IMPORTANT: Material will freeze if stored below 20°C.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000