n-Heptane 95%

Baker

Material No: 8270.9200ST
Effective from: 2013/05/14

Product Specification

Test	Specification
Assay (by GC)	>= 95 %
Boiling Point	98-99°C
Residue after Evaporation	<= 0.001 %

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000