Certified Quality Management System

GOSSELIN™ A Corning Brand

TECHNICAL DATA SHEET

DESCRIPTION : ROUND BOTTLE 100ML HDPE NA H80Ø37 & SCREW CAP RE & SEAL NA, STERILE SAL 10-3

Caractéristiques techniques

- **Type of closure**
  double closure
- **Use temperature of the secondary body**
  -80°C/+120°C
- **Shape**
  assembled
- **Height (mm)**
  80.00
- **Length (mm)**
  51.00
- **Width**
  51.00
- **Diameter**
  37.00
- **Nominal capacity**
  100ML
- **Material and color of the bottle**
  Natural HDPE
- **Moulded graduations**
  N/A
- **Color of the secondary body**
  Red
- **Material of the secondary body**
  HDPE
- **Material of the third body**
  LDPE
- **Physical resistance**
  Flexible, high mechanical resistance
- **Use temperature main body**
  -80°C/+120°C
- **Size label**
  N/A
- **Autoclavage of the main item**
  N/A

Conditionnement

- **Storage Condition (°C)**
  Room température (from 15 to 30°C)
- **Quantity per carton**
  335

Gosselin SAS - 123 Route de Caestre - 59529 Hazebrouck Cedex - France
t +33 328-419-303 - f +33 328-495-692 - communication@gosselin.eu

- **Inner Packaging**
  1 outer bag
- **Packaging**
  by 335 bottles

Références normatives

- **ISO 11137**
  N/A
- **Non-haemolytic**
  N/A
- **Non-pyrogenic**
  N/A
- **Non-cytotoxic**
  N/A
- **Absence of heavy metals**
  Raw material is compliant with the EC/94/62 directive
- **Standard NF en 13972**
  in rigid plastic containers. Nominal capacity definition, brim and total and measurement capabilities brim and total
- **Food contact**
  Raw material is compliant with regulation EC N°1935/2004
- **Based on ISO 11137**
  N/A
- **CE marking**
  N/A
- **Directive 98/79/EC**
  N/A
- **Expiry**
  3 year(s)
- **Sterility Assurance Level based on ISO 11137**
  STERILE SAL 10-3

Précaution d'utilisation

- **Single use**
- **Negative temperature**
  N/A

Domaines d'application

- **Industries**
- **Clinical**
- **Microbiology**

Disponibilité produit

- **Europe Middle East Africa**

- **Asia Pacific**

- **Latin America**

Gosselin SAS - 123 Route de Caestre - 59529 Hazebrouck Cedex - France
t +33 328-419-303 - f +33 328-495-692 - communication@gosselin.eu