Data Sheet

Article: Rührstab 200x5 mm Soda - STIRRING 200x5 Soda Lime Glass

Article no.: 91604005 (L1420005009)

Material: Soda Lime Glass
(Water resistance DIN ISO 719: HGB 3;
Acid resistance DIN 12116: S1;
Alkali resistance ISO 695: A2)

Autoclavable

Dimension/Standard: Length L: 200 mm
Outside diameter D: 5 mm

Handling instructions

The product is made for general use in laboratories.

Due to transport and storage of the products, it cannot be excluded that substances or materials may adhere to the product. Appropriate preparation and cleaning of the products before the planned application is recommended.

Glass will break as a result of impact. Use care when handling to avoid impacting hard objects, such as a spigot, other glassware, counter tops, etc.

Do not evacuate or pressurize unless it is recommended.

Substances can splash out of the glass containers. Wearing safety glasses is recommended.