Potassium Dichromate 1/60 mol/l

DILUT-IT

Material No: 4671
Effective from: 2023/03/30

Product Specification

Test	Specification
Concentration after dilution to 1000 ml in 20°C	c(K2Cr2O7)=1/60mol/l +/-0,2 %

For Laboratory, Research or Manufacturing use.
Volumetric Concentrate, for dilution to 1 l.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000