
REV	DESCRIPTION OF CHANGE	DRAWN BY	DATE	ECN
B	REVISED DRAWING	KW	12/19/11	2011101211112772
C	REVISED LINER MATERIAL	EL	02/05/21	20200121807540742
D	REVISED DRAWING	EL	03/05/21	2021021111121868



FINISH DETAIL
SCALE 2:1

NOTES:
1. DEPRESSED THREADS AT SEAMS.
2. LINER MATERIAL: 1.27 F217 THICK (POLYETHYLENE FOAM) WITH 0.13 THICK PTFE-FACED. LAMINATED TOGETHER. HOT MELT C ADHESIVE TO CAP.
3. BOTTLE CAPACITY: 985.5-10081cc OVERFLOW.
4. WEIGHT: 464.g (BOTTLE ONLY)
5. BOTTOM VIEW MAY VARY DUE TO DIFFERENT MANUFACTURING LOTS.

ITEM	QTY.	MATERIAL: DESCRIPTION AND SPECIFICATIONS
2	12	WHITE POLYPROPYLENE | RESIN: BRASKEM FT-200-WV - SEE NOTE 2
1	12	TYPE III CLEAR SODA-LIME GLASS

THIRD ANGLE PROJECTION

TITLE: 32OZ BTL,BS RD,TYP III AMB,PP CAP,PTFE L

THIS CONFIDENTIAL AND PROPRIETARY DRAWING AND THE INFORMATION IN IT ARE THE EXCLUSIVE PROPERTY OF DWK LIFE SCIENCES. IT IS PROVIDED WITH THE UNDERSTANDING AND INTENTION THAT NOTHING CONTAINED HEREIN WILL BE REPRODUCED OR DISCLOSED TO THIRD PARTIES WITHOUT THE PRIOR EXPRESS WRITTEN CONSENT OF DWK LIFE SCIENCES.

ELECTRONIC OR PAPER COPIES OF THIS PRINT OUTSIDE OF DWK LIFE SCIENCES ARE UNCONTROLLED AND ARE FOR REFERENCE ONLY UNLESS OTHERWISE STAMPED.

TOLERANCES
DIMENSIONS IN MILLIMETERS UNLESS SPECIFIED
X.X ± 2mm
X.XX ± 0.5mm
X.XXX ± 0.13mm

DURAN WHEATON KIMBLE

ITEM NUMBER: W216847	DRAWING NUMBER: CSD-W216847	DRAWN BY: TRJ	DRAWN DATE: 03/28/06		
SAP NUMBER: W216847	OLD MATERIAL NUMBER: N/A	SHEET SIZE: C	SHEET NUMBER: 1 OF 1	SHEET SCALE: 1:1	