1. Part Number: 9286H-100S

2. Photograph:





3. Product Information:

Description:	SteriWare® ULTRA Volumetric Scoop 100ml – Sterile
Use:	Double bagged, single use, disposable powder scoop
Material of Construction:	100% High Density Polyethylene (HDPE) • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • Conforms to USP Class VI • Conforms to USP 661 • Certified BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	100ml
Overall Length	200mm
Handle Length:	95mm
Maximum Width:	59mm
Nominal Weight:	21g
Moulding Environment:	Class 100,000 Cleanroom
Packing Environment:	Class 10,000 Cleanroom

Sampling Systems Ltd.
2 Trillennium Point, Gorsey Lane, Coleshill. B46 1JU. UK.
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

SteriWare is a registered Trademark belonging to Sampling Systems Ltd

Bags:	- Double Bagged, heat-sealed PE/PET bag. - Inner bag printed with batch information
Outer Packaging:	Cardboard box
Box Liner:	LDPE bag
Quantity per Box:	50
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to ***change.control@sampling.com***

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable product. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 8-Jan-25

Sampling Systems Ltd.
2 Trillennium Point, Gorsey Lane, Coleshill. B46 1JU. UK.
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

SteriWare is a registered Trademark belonging to Sampling Systems Ltd