B BRAUN
SHARING EXPERTISE

PRODUCT DATA SHEET

PERFUSOR LINE, PE, LL, 100 CM (8255067)

For use with syringe pumps



Attribute Name	Attribute Value
Product type	Administration set
Brand name	Perfusor®
Product name (HC)	Original Perfusor® Line
Type of filter	No filter
Position of filter	-
Number of SafSite®	-
Number of needle-free valves	-
Number of back check valves	-
Tube length	100.00 cm
Outer-ø tube	2.00 mm
Inner-ø tube	1.00 mm
Priming volume	0.880 ml

B. Braun Melsungen AG
Carl-Braun-Straße 1
34212 Melsungen | Germany
Tel. +49 5661 710
www.bbraun.com

Attribute Name	Attribute Value
max. pressure	2.0 bar
Connection	Luer-Lock/male, Luer-Lock/female
Light-protected	No
Color of light protection	-
with Cannula	no
Gauge (cannula)	-
Material	PE
Type of packaging	Blister
Product (unpack.) contains PVC	No
Product (unpack.) cont. Latex	No
Product (unpack.) cont. DEHP	No
Packaging contains PVC	No
Packaging contains Latex	No
Packaging contains DEHP	No
Product contains PVC	No
Product contains Latex	No
Product contains DEHP	No
Product contains Nickel	No
Sterility	sterile
Type of Sterilization	EO
Type of Sterilizat. (allowed)	EO
Disposable / Reusable	disposable

B. Braun Melsungen AG
Carl-Braun-Straße 1
34212 Melsungen | Germany
Tel. +49 5661 710
www.bbraun.com