LP ITALIANA SPA
Via C. Reale 15/4, 20157 Milano, ITALY
T +39 02 3933 06.1 - F +39 02 3931 3484
www.lpitaliana.com - info@lpitaliana.com

Capitale Sociale € 309.600,00
R.E.A. MI 882798
Reg. Imp. C.F. e P.I. MI 01794050151

PRODUCT DATA SHEET

SPECIMEN TRANSPORT BAG

LP ITALIANA's low density polyethylene double-pouch bags are designed for transporting a patient's biological samples together with the relative documents to the laboratory.

Code	**Colour**	Dimensions (cm)	
		Pouch with closure	Document pouch
172030	Orange	15 x 17	15 x 23
172031	Red		
172033	Yellow		
172034	Green		
172035	Light blue		
172037	Neutral		
172045	Pink		
172130	Orange	16 x 20	16 x 28
179137	Neutral	20,5 x 32,5	20,5 x 38
179237	Neutral	30 x 42	//

Main features:

1) The pouch with mini-grip type closure ensures a leak-proof closure, in case there is some sample leakage from the enclosed test tube(s) or container.
2) The document pouch and sealable pouch are permanently attached, to ensure that samples and documents are kept together during transport, thus avoiding mix-ups.
3) LP transport bags may also be used in pneumatic delivery systems, due to the closure.
4) An orange coloured Biohazard symbol is printed on the bags.
5) A white writing area is located on the bags.
6) Bags are available in various colours for easy identification and sorting.

All dimensions indicated are nominal. Colours may vary slightly between production lots. Technical specifications are subject to change without prior notice.