Technical Datasheet

VWR WIPES CLEANSMOOTH TURQUOISE

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Applicable part numbers : 115-2396, 115-2397

Attribute	Value*	Test method
Composition :		
wood pulp	55%	
polyester	45%	
Production method	spunlace	
Colour	turquoise	
Basis weight	73 g/m²	
Absorption	600%	
Tensile strength :		
Machine direction	146,8 N	
Cross direction	80,1 N	

Characteristics :
Low fiber and particle release ; Very strong, even wet ; Resistant to acids and solvents ; Good absorption. Recyclable. Food Contact Certified. Does not contain binders. Silicone-free.

Applications :
Pharma ; Aerospace ; Semiconductor ; Food & Beverage, Automotive, Laboratories ; Printing, Healthcare.

Formats & Packaging :

Cat No	Type	Format size	Quantity	Primary pkg	Secondary pkg	Pallet
115-2396	Industrial roll	38x30cm	2 rolls x 500 pcs	Plastic film	n/a	42 x 2rolls
115-2397	Z-folded wipes	38x30cm	10 packs x 50 pcs	Plastic bag	Cardboard case	64 cases

**** Disclaimer:*** *The information presented here is the result of testing and other evaluations conducted during the development of the product. This information does not constitute a specification.*

Date of issue: **March 2017**