You are on MN's World site

Stay on World site

Do you want to choose another region?

EN

Hardware	Polypropylene (PP) housing with filter membrane
Membrane type	Polytetrafluoroethylene (PTFE)
Membrane properties	Hydrophobic
Membrane diameter	25 mm
Membrane pore size	1 µm
Code type	Labeled
Color code (top)	Labeled
Color code (bottom)	Colorless
Inlet	Luer lock (top)
Outlet	Luer tip (bottom)
Sample volume	5-100 mL
Recommended application(s)	Nonpolar samples, Suitable for HPLC columns packed with particles > 3 µm, GC, SFC, ..., Very resistant towards all kinds of solvents as well as acids and bases
Sterile	No
Binding capacity (proteins)	Moderate–high
Dead volume	80 µL
Filtration area	4.91 cm^2
Storage temperature	RT

Related products +