VWR We Enable Science

TECHNICAL DATA SHEET


PLASTIC CANISTER 5 LITER (MARWIN Type)


SPECIFICATIONS

Product:	Stackable Canister 5 Liter
Material:	PE-HD
UN-certification:	3H1/Y1.9/200/../A/
Closure:	DIN 51
Dimensions (LxWxH):	248 x 142 x 217 mm









Remark : Packaging materials and specifications shown may be modified without prior notice
Version dd. June 24, 2020

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343