B BRAUN
SHARING EXPERTISE

PRODUCT DATA SHEET

DISCOFIX C MANIFOLD 5 (MULTI) (16608C)



Attribute Name	Attribute Value
Brand name	Discofix®
Product name (HC)	Discofix® C
Product type	Manifold
Type of stopcock manifold	5-fold
Accessories	No
Color	multicolored
Lipid resistant	Yes
Drug resistant	Yes
max. pressure	2.0 bar
Type of packaging	Peel pack
Product (unpack.) contains PVC	No
Product (unpack.) cont. Latex	No
Product (unpack.) cont. DEHP	No

B. Braun Melsungen AG
Carl-Braun-Straße 1
34212 Melsungen | Germany
Tel. +49 5661 710
www.bbraun.com

Attribute Name	Attribute Value
Packaging contains PVC	No
Packaging contains Latex	No
Packaging contains DEHP	No
Product contains PVC	No
Product contains Latex	No
Product contains DEHP	No
Sterility	sterile
Type of Sterilization	Gamma
Disposable / Reusable	disposable

B. Braun Melsungen AG
Carl-Braun-Straße 1
34212 Melsungen | Germany
Tel. +49 5661 710
www.bbraun.com