# DURAN® **YOU**TILITY IDENTIFICATION **LABELS**





The accurate and reliable labelling of scientific samples is an integral part of the success of your laboratory. Labels carry critical information, and without accurate labelling, months of research could be put at risk. Identification labels are the foundation of an effective laboratory management process.

The robust DURAN® **YOU**TILITY identification labels are suitable for use with the **YOU**TILITY bottle, plus many other types of laboratory glassware, such as GL 45 bottles, beakers and conical flasks.

DURAN
WHEATON
KIMBLE

# YOUR DESIGN:



DESIGN




PRINT



STICK

Create and print out your own custom label designs online with our new DURAN® label creator web app:

## DURANLABELS.COM

- The tear resistant, white polyester **YOU**TILITY labels use a high performance, peelable adhesive.
- An easy-to-use label creator web app. is available (www.duranlabels.com) to design, and print your own label designs.
- Ideal for use in fridges, freezers, autoclaves, incubators, and water baths, without the risk of the labels falling off.
- Can be printed with office printers, copiers, or marked with technical lab pens (pre-test ink for suitability).
- Chemically resistant to many typical laboratory chemicals, and solvents.
- Wide –40 to +150 °C thermal performance range.
- Supplied in a re-sealable, pack of 100 labels on five A4 sheets.

| Cat. No. | | DESCRIPTION | MATERIAL AND COLOUR | PACK SIZE |
|---|---|---|---|---|
| NEW | 29 401 02 03 | DURAN® YOUTILITY Printable Identification labels | White polyester 36 x 70 mm | 1 x 100 labels |

You can obtain this product from your specialist labware distributor. Any product related or technical questions? Get in touch with us by e-mail: sales@duran-group.com

DURAN
WHEATON
KIMBLE

**DWK Life Sciences GmbH**
Otto-Schott-Str. 21
97877 Wertheim / Main
Germany

Phone: +49 (0) 6131 1445 4131
Fax: +49 (0) 6131 1445 4016
sales @ DWK-LifeSciences.com
www.DWK-LifeSciences.com