N,N-Dimethylformamide

AR

Material No: 6785-25
Effective from: 2013/03/25

# Product Specification

| Test | Specification |
|---|---|
| Assay (GC) | >= 99.8 % |
| Water (KF) | <= 0.07 % |
| Titrable Acid (meq/g) | <= 0.0005 |
| Residue after Evaporation | <= 0.005 % |

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000