DATA SHEET Gr. 0606

●DESIGNATION

-Bench protective paper
-Material: Cellulose + Polyethylene film
-Binder: no

●APPLICATIONS

-Protection of laboratory benches, working areas cupboard or trays
-Protection against contamination, chemical pills, toxic infectious and radioactive materials
-Suitable for laying out of animal cages

●TECHNICAL SPECIFICATIONS

Characteristics	Values	Unities	Standards
Weight	225 (+/- 5%)	g/sqm	150536
Thickness	0,50	mm	
Capillary rise	88-93	mm	Klemm
Capacity of absorption	400	g/m²	