boettger

Technical Data (TD)

Tube PS 12/75 mm, Round					
REF	01-051-0100				

1.	**Description / Specification**			
1.01	Description	Tube 12/75 with round bottom, appropriate grip stopper REF 01-624-010x		
1.02	Dimension	See *Customer Drawing*		
1.03	Volume	5 ml		
1.04	Material / Resin	Polystyrene (PS)		
	Grade of Resin	MAT100		
	Properties of Resin	See *www.boettger.com*		
1.05	Colour	Crystal Clear		
1.06	Hygiene Level	Standard hygienic		
1.07	Quality Control	Testing of attributive and variable characteristics in accordance with the valid specification		
1.08	Storage	To be stored in closed original packaging: Dry, at room temperature		
1.09	Usage / Purpose	General laboratory product, Collection and evaluation of samples		
1.10	Conformities / Regulations	See *Classification*		
1.11	Product Stewardship	See *Product Stewardship* of Grades of Resins/Pigments		
1.12	Other information	For single use only		
2.	**Features**			
2.01	Basic features	-		
2.02	Temperature range	-10°C to 60°C		
2.03	Autoclavability	No		
2.04	Centrifugation	2000 g		
2.05	Chemical Resistance	See *www.boettger.com*		
2.06	Expiry	5 years from LOT (YYYY-MM)		
2.07	Other information	-		
3.	**Packaging**			
3.01	Pcs/Bag	1000		
3.02	Pcs/Ctn	4000		
3.03	Pcs/Pallet	80000		
3.04	LOT / VP	WWYY (week, year of production) / XX (packaging)		
3.05	Other information	See *Measures*		
4.	**Other Information**			
4.01	Codes and Numbers	See *Codes and Numbers*		
4.02	Symbols	See *Symbols*		

All informations without guarantee. Above mentioned products are manufactured according to company specifications and our certified Quality Management System corresponding to DIN EN ISO 13485. The quality check is carried out on a sample basis; an individual examination is not possible. Deficiencies parts may occur. We accept no liability for consequential damages.