TECHNICAL INFORMATION SHEET

BD Vacutainer® Urine Transfer Straw

Product Catalogue Number: 364940

Product Description

The Urine Transfer Straw is a single use, non-sterile device intended for the transfer of urine from a primary containment vessel to a urine collection tube. BD claims that use of this device eliminates the need for pouring. This product is intended for use by healthcare professionals

Manufacturing Information

(Legal) Manufacturer:	Becton Dickinson and Company, 1 Becton Drive, Franklin Lakes, NJ 07417, USA
Standards & Certificate Numbers:	ISO 13485, MD 19.2137
Country of origin:	USA
Notified Body:	NSAI (0050)
Certification body:	NSAI (0050)
EU Authorised Representative:	BD Switzerland Sarl, Terre Bonne Park - A4, Route de Crassier 17, 1262 Eysins, Switzerland

Sterilisation

Method:	This product is non-sterile

Compliance

Directive:	European In Vitro Diagnostic Medical Devices Directive 98/79/EC
Classification:	Non Annex II

Product Specification

Product Storage:	Do not expose to direct sunlight
Shelf-life:	No specified shelf life
Cannula Material:	Stainless Steel (304 Grade)
Holder Material:	High Density Polyethylene (HDPE)
Straw Material:	Polystyrene
Cannula Lenght (mm):	39.3 mm
Latex (NRL):	No
Dry Natural Rubber (DNR):	No
Phthalates:	No
Material of animal origin:	No



Packaging Specifications

100 unit pack weight (kg):	0.3497	100 unit packaging material:	Plastic Bag
100 unit pack volume (m^3):	Unknown	100 unit packaging weight (kg):	0.016667
100 unit pack dimensions LxHxW (mm):	320 x 130 x 255	1000 unit pack weight (kg):	3.497
1000 unit packaging material:	Cardboard	1000 unit pack volume (m^3):	0.0504
1000 unit packaging weight (kg):	0.657	1000 unit pack dimensions LxHxW (mm):	400 x 450 x 279
1000 unit packaging volume (m^3):	Unknown		

Product Catalogue Number: **364940**

Labelling Information

All labelling complies with the requirements of the European In Vitro Diagnostic Medical Devices Directive 98/79/EC and includes the CE marking.

	Unit Pack	Shelf Pack	Case Pack
Company name		•	•
Reference number		•	•
CE marking		•	•
Single use symbols		•	•
Lot number		•	•
Primary barcode (GS1-128) product identification		•	
Storage instructions		•	•
Quantity in package		•	•
Product name & short description		•	•
Authorised Representative		•	•

Instructions For Use



Further Reading

Not applicable

Sample Storage & Stability

Not applicable

References

Not applicable

© 2020 BD, BD logo and all other trademarks are the property of Becton, Dickinson and Company. 364940