
LE407

pH ● Ref ● °C ○ ISM

For general applications;
robust and maintenance-free



Sensor type	combined pH electrode
Order number	51340330
Specifications	
pH range	0...14
zero point (pH)	7.00 ± 0.25
temperature range	0...80 °C
temperature probe	-
shaft material	POM
membrane glass	U
membrane resistance	< 250MΩ
reference system	Ag / AgCl
reference diaphragm	Polyester
reference electrolyte	gel
bridge electrolyte	-
electrode head	-
cable length	1.0 m
connector	BNC
cell constant (cm^{-1})	-
cell type	-
indicating metal	-
Storage	
3M KCl (51350072)	

© 01/2018 MTII

METTLER TOLEDO