Cole-Parmer®
stuart™

Cole-Parmer® Stuart™ Incubator – INC-200 & INC-200D

- Incubation of complete instruments
- Choice of analogue or digital temperature control
- Easy access to working chamber
- Full visibility of samples
- Suitable for multiple applications



Incubator INC-200

coleparmer.com

antylia SCIENTIFIC



Incubator INC-200D

Cole-Parmer® Stuart™ Incubator – INC-200 & INC-200D

The Cole-Parmer® Stuart™ INC-200 & INC-200D Incubators are suitable for multiple applications such as plant propagation, humidity testing, climactic condition simulation, corrosion testing and more.

Forced air circulation and electronic control ensures accurate temperature conditions. A pre-set push switch automatically sets the temperature at 37°C, while a safety feature cuts off heating if the unit reaches 72°C or higher.

Temperature is set via a calibrated knob for the analogue model (INC-200), and via an LCD display when using the digital model (INC-200D). The latter not only allows more accurate temperature management but also gives a constant reading of actual temperature when in use.

The purpose-built design offers easy access to samples - a hinged front door panel and plugged holes (2 x 10mm diameter) on each of the side panels for the introduction of gases or cables. The clear acrylic® polymer construction also allows total sample visibility at all times.

Technical Specification

Specification	INC-200/INC-200D
Temperature range	Ambient +5°C to 60°C
Temperature fluctuation at 37oC	±0.1°C
Temp. variation between shelves	±0.3°C
Nominal volume	60 litres
Internal dimensions, mm (w x d x h)	450 x 380 x 380
Overall dimensions, mm (w x d x h)	600 x 390 x 390
Net weight, kg	11.2
Electrical supply	230V, 50Hz, 350W
IP Rating	30

Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Acrylic incubator analogue	01110-02	INC-200	INC-200	SI60
Acrylic incubator analogue 120V/60HZ	01110-00	INC-200	INC-200-120	SI60/120V/60
Acrylic incubator digital	01110-08	INC-200	INC-200D	SI60D
Acrylic incubator digital 120V/60HZ	01110-04	INC-200	INC-200D-120	SI60D/120V/60

Accessories

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Acrylic base plate	01110-97	INC-200	N/A	5160/1
Plastic coated shelf/rack system with two shelves	01110-98	INC-200	N/A	S160/2

USA:
+1.800.323.4340
+1.847.549.7600

Canada: +1.800.363.5900
China: +86.21.5109.9909
France: +33 (0) 1486 37800
Germany: +49 (0) 9377 92030

India: +1.800.266.1244
Italy: +39 (0) 1313 89513
UK: +44 (0) 1480 272279
All other countries: +1.847.549.7600

Cole-Parmer® essentials
an antylia scientific company

Find out more at antylia.com

AZ-EMEA-LT-0018 0721