3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Bromoform - contains 60-120 ppm 2-methyl-2-butene as stabilizer, 99%

Product Number:	**241032**
CAS Number:	75-25-2
MDL:	MFCD00000128
Formula:	CHBr3
Formula Weight:	252.73 g/mol



TEST	Specification
Appearance (Color)	Colorless to Light Yellow
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
GC (area %)	≥ 98.5 %
Stabilizer	Confirmed
Contains 60-120 ppm 2-Methyl-2-butene	

Specification: PRD.1.ZQ5.10000003769

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.