Chloroform

AR

Stabilized with Amylene

Material No: 5411-45
Effective from: 2014/03/17

Product Specification

Test	Specification
Assay (GC)	>= 99.5 %
Water (KF)	<= 0.05 %
Titrable Acid (meq/g)	<= 0.0005
Residue after Evaporation	<= 0.001 %
Stabiliser (Amylene)	~50 ppm
Free Chlorine (Cl2)	<= 0.0005 %

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000