1. Product name:

Teflon connector.

2. Catalogue No.:

1 71838.

3. Drawing No.:

3033

4. Product application:

The product is used as a connector piece between cell factories cover caps and tubes. The product is mounted with two O-rings, in the end that is put down into the cover caps.

5. Description:

5.1 Material:

Teflon, White. Recycling code 7 according to SPI-norm.

5.2 Dimensions and tolerance:

See drawing 3033

5.3 O-rings:

Material: Silicone

5.4 Product property:

Autoclavable

6. Delivery:

6.1 Packing and packaging:

10 connectors are put in a double bag.

6.2 Labelling:

The bag are labelled with product No., Lot No. and number of units.