SciLabware Limited
Unit 4, Riverside 2, Campbell Road
Stoke-on-Trent, ST4 4RJ, United Kingdom
Tel +44 (0) 1782 444406
Fax +44 (0) 1782 940436
www.scilabware.com
www.wheaton-uk.com
Registered office (as above)
Registered in England and Wales No. 7723574
BS EN ISO 9001: 2008 registered company

QF/7.17 Issue A 10 Nov 2017

PRODUCT SPECIFICATION

Product Code: 1003-1L

Product Description: 1000ml Heavy Duty, Low form Beaker

Brand: PYREX

Type of Glass: Borosilicate glass Type 3.3, which conforms to ISO3585 European and US Pharmacopoeia Type 1 glass

Product Specification: White graduation and marking spot

Height	156mm
O.D.	108mm

Pack Quantity: 6 pieces per pack

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001 and BS EN ISO13485.

Signed on behalf of SciLabware Limited by:

Signed: M.Ashley
Quality Technician

Date: 22 Oct 2019

Pyrex® is a registered trademark. SciLabware registered user.