
Ed.: 1 . Vig.: 15.01.2011 . Prod.: 283462

PanReac AppliChem ITW Reagents

PRODUCT CODE: 283462

Ferroin solution 0.025 mol/l (0.025M) for volumetric analysis
redox indicator

$C_{36}H_{24}FeN_6O_4S$
$C_{36}H_{24}FeN_6O_4S$



M.= 692,24 CAS [14634-91-4] EINECS 238-676-6
TARIC 3822 00 00 00

SYNONYMS: 1,10-Phenanthroline-Iron(II) Sulphate Complex, Tris (1,10-Phenanthroline) Iron(II) Sulphate

PHYSICAL DATA: Miscible with water D 20/4 1,008 • pH(50 g/l)2,5 •
BIBLIOGRAPHY: Merck Index **12,** 7.356

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l~1,008 kg 1kg~0,992 l

SPECIFICATIONS:

Suitability: as redox indicator	passes test

Ed.: 1 . Vig.: 15.01.2011 . Prod.: 283462