Acetonitrile

Baker Ultra Resi-Analyzed

Material No: 9255-02
Effective from: 2017/04/04

Product Specification

Test	Specification
Assay (by GC) (corrected for water)	>= 99.8 %
Color (APHA)	<= 10
Residue after Evaporation	<= 1.0 ppm
Titrable Acid (ueg/g)	<= 0.3
Titrable Base (ueg/g)	<= 0.6
Water (H2O)	<= 0.05 %
ECD Sensitive Impurities (as Heptachlor Epoxide)	
Single Impurity Peak (pg/ml)	<= 10
FID-Sensitive Impurities (as 2-Octanol)	
Single Impurity Peak (ng/ml)	<= 5

For Laboratory, Research or Manufacturing use.
For Organic Residue Analysis.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000