Technical Data Sheet



VWR® 500ml Centrifuge Bottles

Description: 500mL, Conical Centrifuge Bottle, Flat Cap, Sterilized


Purpose: Centrifuge bottles are popular in samples centrifuging or storage in many research area

Materials:

Bottle: PP(Polypropylene), Compliance with USP Class VI

Cap: HDPE(High Density Polyethylene), Compliance with USP Class VI

Dimensions (Unit: mm)



Features

- Made from high grade gamma resistance polypropylene
- Flat cap made with polyethylene
- Meet the requirements of BSE/TSE
- Latex Statement: this product is latex free
- Engraved graduation
- Max RCF up to 6,000xg
- D/RNaseNase free and non-pyrogenic
- Sterilized by irradiation SAL10^{-6} (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system

Centrifuge Adapter

No.	Brand Name	Type	Rotors
1	BECKMAN	J-E	JLA-9.1000

② STERILE | R

VWR® 500mL Centrifuge Bottles

VWR EU EAN. No.	VWR NA Cat. No.	Capacity (mL)	Maximum RCF	Bottle Color	Cap Color	Cap Shape	Bottle Shape	Sterile	Packaging Configuration
525-1598	76234-746	500mL	6000g	Clear	White	Flat cap	Conical-Bottom	Yes	6/Bag, 36/CS