Product Name: Potassium chloride, 0.010M, Conductivity Standard, Thermo Scientific Chemicals

Catalog Number: 042679.AP

CAS Number: 7447-40-7

Molecular Formula: ClK

Molecular Weight: 74.55

MDL Number: MFCD00011360

InChI Key: WCUXLLCKKVVCTQ-UHFFFAOYSA-M

SMILES: [Cl-].[K+]

Synonym: potassium chlorideIsuper KImuriate of potashIklor-conIsylviteICIKIKay-EMIKCIImicro-KIslow-K

Product Specification

Conductivity: 1413micromho at 25?C

Concentration: 0.010M

Date Of Print: 10/29/2023

Product Specifications are subject to amendment and may change over time. Data contained is accurate as of the date printed.