Certificate

Filter paper MN 280

REF: 100L280001A

Product description

ashless, slow filtration

typical data	unit	specifies	based on
basis weight	g/m²	75	DIN EN ISO 536
thickness	mm	0.15	DIN EN ISO 534
filtration speed	s	100	DIN 53137
ash content	%	< 0.01	DIN 54370
surface		smooth	

Confirmation

Hereby we confirm that the above mentioned product has successfully passed our quality control system in accordance with EN ISO 9001 and meets the specific quality criteria.

This document has been produced electronically and is valid without a signature.

MACHEREY-NAGEL GmbH & Co. KG
Valencienner Str. 11
52355 Düren · Germany
www.mn-net.com

DE Tel.: +49 24 21 969-0 info@mn-net.com
CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com
FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com
US Tel.: +1 888 321 62 24 sales-us@mn-net.com