20 µl Bevelled TipOne® Filter Tip, Natural, Racks (sterile)

Art. No. S1120-1810



Product Highlights

- Highly effective HDPE filter that contains no additive to interfere with your samples
- Universal fit. One tip fits all! The tip collar is
- engineered to fit all popular pipettes.
- Highly polished moulds for an extremely smooth surface to reduce sample retention
- Highest quality product. Tips manufactured from pure, virgin polypropylene to ensure the highest quality product.
- Certified RNase, DNase, DNA and Pyrogen free.

Product Applications

DNA applications (e.g. PCR)
RNA applications (e.g. Gene expression analysis)
Protein applications (e.g. Antibody Research)
Cell Culture applications (e.g. Media)
Applications with radio actives

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

STARLAB (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk

General Data

Pack Size	Rack
For use with	Abimed, Biohit®, BRAND®, Capp®, Eppendorf®, Gilson®, HTL, Nichiryo®, Rainin®, Sartorius®, Socorex®, STARLAB®, Thermo Scientific™
Color Code Wafer	Yellow
Minimum volume	1.0 µl
Material	Polypropylene
Color	Natural
Autoclavable	No
Sterile	Yes
Volume	20 µl
Packaging Item	Tips
Volume range	1.0 – 20 µl
Filter tip	Filter
Packaging Unit	Racks
Pack Size	960 Tips (10 Racks × 96 Tips)
For use with	Abimed, Biohit®, BRAND®, Capp®, Eppendorf®, Gilson®, HTL, Nichiryo®, Rainin®, Sartorius®, Socorex®, STARLAB®, Thermo Scientific™

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

STARLAB (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk

For use with	**Abimed:** BioPette-E, Discovery Comfort, LABMATE Optima **Biohit®:** eLine, mLine, Picus, Proline **BRAND®:** Transferpette, Transferpette Electronic, Transferpette S **Capp®:** CappAero Comfort, CappAero Universal, Capptronic Multi, Capptronic Single, ecopipette, Singlechannel **Eppendorf®:** Reference, Reference 2, Research, Research plus **Gilson®:** Pipetman Neo, Pipetman P, Pipetman Ultra **HTL:** Discovery Comfort, Discovery Pro, Labmate Pro, Optipette **Nichiryo®:** 50F, NLE, NPLO, NPX **Rainin®:** Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus **Sartorius®:** eLine, mLine, Picus, Picus NxT, Proline, Proline Plus, Tacta **Socorex®:** Acura 825, Calibra digital **STARLAB®:** ErgoOne, ErgoOne E, MicroOne, StarPet **Thermo Fisher Scientific Inc.:** Finn Focus, Finn Focus Univ., Finnpipette Digital, Finnpipette F1/F2, Finnpipette Novus, Finnpipette Novus Yellow,

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

STARLAB (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk

More informations about 20 µl Bevelled TipOne® Filter Tip, Natural, Racks (sterile)

Nothing goes through!

Protect your samples against unwanted contamination by using TipOne® filter tips. The TipOne® filter is one of the best on the market and protects against aerosol cross-contamination without trapping your valuable samples. The filters contain no additives to interfere with samples.

During pipetting it is not possible to see with the eye if fine aerosols are going into your pipette. If they are, those aerosols can contaminate other samples at a later date.

TipOne® filter tips are proven to block >99% of aerosols and are among the best filters on the market. The results, shown right, speak for themselves!

**The comparative tests were conducted in an independent test institute with filter tips from other leading manufacturers.*

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

STARLAB (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk