Seit 1927
über 90 Jahre

Since 1927
more than 90 Years

Amarell
GmbH & Co. KG
Präzisions-Thermometer / Aräometer

Lindenstraße 43
D-97892 Kreuzwertheim

Tel. +49 (0) 9342 9283-0
Fax: +49 (0) 9342 39860

Internet: www.amarell.de
E-Mail: thermometer@amarell.de

Technical Data Sheet

Art. No.:	L33004-TOL
Type:	Precision thermometer, low temperature, DIN 12778, enclosed scale, cap. lens front, colourless
Measurement range:	-50+50°C
Subdivision:	1°C
Length:	300mm
Diameter:	8-9mm
Filling:	Toluene
Immersion:	Total
Finish:	Richter-top finish with button
Certificate:	Without Certificate

August'23