SPECIFICATION OF 12" POWDER FREE VINYL GLOVES

(AMBIDEXTROUS TYPE, NATURAL COLOR, SMOOTH TYPE)

SIZE / SPECIFICATION	X-SMALL	SMALL	MEDIUM	LARGE	X-LARGE	TEST CODE
WIDTH OF PALM (m/m)		88 ± 2	98 ± 2	108 ± 2	113 ± 2	ASTM D5250
TOTAL LENGTH (m/m) from middle finger tip to cuff	min. 300					
LENGTH OF FINGERS(m/m)						± 1.0
THUMB		56	61	62	67	
INDEX FINGER		64	72	77	85	
MIDDLE FINGER		73	83	91	97	
RING FINGER		69	77	83	89	
LITTLE FINGER		53	59	62	66	
CIRCUMFERENCE OF FINGERS (m/m)						± 1.0
THUMB		66	74	82	82	
INDEX FINGER		58	62	70	71	
MIDDLE FINGER		56	60	70	71	
RING FINGER		54	60	70	70	
LITTLE FINGER		50	56	63	63	
Before and after accelerated aging (70 ±2℃ x 72 ±2h) TENSILE STRENGTH(MPa)	min. 11 Mpa					ASTM D5250
Before and after accelerated aging ELONGATION (%)	min. 300 %					ASTM D5250
THICKNESS (m/m)	FINGER TIP min. 0.05					
	PALM min. 0.08					