Wide-necked wash bottle



Wash bottle with medium overprint as well as hazard and safety instructions (health statements/physical statements) based on the GHS system with substance name in German, English, French and Spanish, GHS pictograms and CAS number.


- Safe filling due to wide-necked opening 39 mm (neck dia. inside)
- LDPE yellow transparent with light protection
- Contents 500 ml

LaboPlast

Wide-necked wash bottle			
Print	Colour	Item no.	
without label - neutral	yellow	0310-3050	N3
Distilled water	transparent	0310-3051	N3
Methanol	green	0310-3052	N3
Ethanol	orange	0310-3053	N3
Isopropanol	blue	0310-3054	N3
Acetone	red	0310-3055	N3

Safety wash bottle





The vapour-permeable membrane reliably prevents leaking of the contents.

If the room temperature increases, or the sun shines into the lab, conventional wash bottles will drip because of the increased inner pressure. The venting valve ensures that the excess pressure can escape and dripping is prevented. With leakproof venting valve, contents cannot leak when bottle has fallen over as with wash bottles with ball valve or venting hole.

- Safe filling due to wide-necked opening 39 mm (neck dia. inside)
- Permanent print according to GHS with substance name in German, English, French and Spanish, GHS pictograms and CAS number
- LDPE yellow transparent with light protection
- Contents 500 ml

LaboPlast

Safety wash bottle			
Print	Colour	Item no.	
without label - neutral	yellow	0310-2056	N3
Distilled water	transparent	0310-2050	N3
Methanol	green	0310-2052	N3
Ethanol	orange	0310-2054	N3
Isopropanol	blue	0310-2055	N3
Acetone	red	0310-2051	N3