all-glass syringes

Well suited for viscous liquids, corrosive reagents, dyes, etc.

Dosys™ *all-glass* 155

The reusable borosilicate glass syringes stand for superior chemical, heat and shock resistance in various application fields.

A precious, low cost alternative to disposable plastic syringes for easy infusing, diluting, aspirating and sampling (not gas tight).



160°C

** New **
JUMBO SIZES
up to 250 mL

Product advantages

- Precision-machined plunger and barrel
- Excellent fitting and liquid tightness
- Permanent, high visibility graduations
- Glass Luer and metal Luer Lock nozzles
- Sterilisation up to 160°C / 320°F
- Very easy to disassemble and clean

Ordering information - Syringes

Volume	Division	Packaging	Glass Luer nozzle Cat. No.	Metal Luer Lock nozzle, Cat. No.
0.1 - 1 mL	0.05 mL	3 / pk	**155.0301**	**155.0501**
0.5 - 2 mL	0.1 mL	3 / pk		**155.0502**
0.2 - 5 mL	0.2 mL	3 / pk	**155.0305**	**155.0505**
1 - 10 mL	0.2 mL	3 / pk	**155.0310**	**155.0510**
1 - 20 mL	1 mL	2 / pk	**155.0320**	**155.0520**
1 - 30 mL	2 mL	2 / pk	**155.0330**	**155.0530**
1 - 50 mL	2 mL	1 / pk		**155.0550**
10 - 100 mL	10 mL	1 / pk	**155.03100**	**155.05100**
10 - 150 mL	10 mL	1 / pk		**155.05150**
10 - 200 mL	10 mL	1 / pk		**155.05200**
10 - 250 mL	10 mL	1 / pk		**155.05250**

Dispensing cannulas

Blunt end, superior stainless steel, Luer Lock, autocl.

Description (Ø x l)	Packaging	Cat. No.
0.8 x 40 mm	12 / pk	**370.0840**
1.2 x 50 mm	1 / pk	**370.1250**
2.2 x 100 mm	1 / pk	**370.22100**
2.2 x 150 mm	1 / pk	**370.22150**

Barrel cleaning brushes

Heavy duty nylon, tufted end and plastic handle

Brushes for syringes	Packaging	Cat. No.
1 mL	6 / pk	**1.170.301**
2, 5 mL	6 / pk	**1.170.305**
From 10 to 250 mL	3 / pk	**1.170.320**
Set of two each	6 / pk	**1.170.300**

For information on our large choice of needles and drenchers or for an overview of the complete Socorex product range, visit www.socorex.com

XF. 155.E – C1012

SOCOREX SWISS
SOCOREX ISBA S.A. ▫ CHAMP COLOMB 7 ▫ P.O. BOX ▫ 1024 ECUBLENS ▫ SWITZERLAND
E-MAIL socorex@socorex.com ▫ www.socorex.com ▫ TEL. +41 21 651 6000 ▫ FAX +41 21 651 6001

G.M. SYSTEM ISO 9001/13485 CERTIFIED