Tetrahydrofuran
ChromAR HPLC

Material No: 2858-25
Effective from: 2013/03/08

# Product Specification

| Test | Specification |
|---|---|
| Assay (GC, corrected for water) | >= 99.8 % |
| Water (KF) | <= 0.03 % |
| Acidity | <= 0.0005 |
| Residue after Evaporation | <= 0.0005 % |
| Peroxides (as H2O2) | <= 0.05 % |
| Ultraviolet Absorbance (1.00-cm cell vs. water): | |
| 240 nm | <= 0.3 |
| 250 nm | <= 0.15 |
| 260 nm | <= 0.08 |
| 270 nm | <= 0.03 |
| 280 nm | <= 0.01 |
| Fluorescence (as quinine) | |
| 254 nm | <= 1 ppb |
| 365 nm | <= 1 ppb |

For Laboratory, Research or Manufacturing use.
Suitable for Liquid Chromatography and UV-Spectrophotometry

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000