Tetrahydrofuran (Stabilized)
HPLC
For Use in Gel Permeation Chromatography
Contains Butylated Hydroxytoluene as a Preservative

Material No.: 9440-03
Revision No.: 0

Product Specification

Test	Specification
Assay (by GC, corrected for water, exclusive of preservative)	≥ 99.5 %
Color (APHA)	≤ 10
Peroxide (as H_2O_2)	≤ 0.005 %
Stabilizer (BHT)	75 - 400 ppm
Water (by KF, coulometric)	≤ 0.010 %

For Laboratory,Research,or Manufacturing Use
Stabilized with approximately 0.01% Butylated Hydroxytoluene

Packaging Site: Paris Mfg Ctr & DC

James T Ethier

Jamie Ethier
Vice President Global Quality

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700
Avantor Performance Materials, LLC
100 Mansford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700