Ethylene Glycol
BAKER ANALYZED® Reagent
(1,2-ethanediol)

Material No.: 9300-03
Revision No.: 0

Product Specification

Test	Specification
Appearance	Passes Test
Assay ($HOCH_2CH_2OH$) (by GC)	≥ 99 %
Assay ($HOCH_2CH_2OH$)(wet analysis)	≥ 99 %
Color (APHA)	≤ 10
Specific Gravity at 20°/20°C	1.114 - 1.116
Residue after Ignition	≤ 0.005 %
Acidity (as CH_3COOH)	≤ 0.01 %
Water (H_2O)(by Karl Fischer titrn)	≤ 0.2 %
Chloride (Cl)	≤ 5 ppm
Trace Impurities - Iron (Fe)	≤ 0.2 ppm

For Laboratory,Research,or Manufacturing Use

Packaging Site: Paris Mfg Ctr & DC