

NOTES:

1. APPROX. WEIGHTS: BOTTLE: 3.0 OZ. | CAP: 4.2 GRAMS AVG EXCLUDING LINER.
2. DEPRESSED THREAD AT SEAMS.
3. SIGHT "H" DOWN "E".
4. DIMENSIONS NOT TOLERANCED ARE FOR REFERENCE ONLY.
5. BOTTOM VIEW MAY VARY DUE TO DIFFERENT MANUFACTURING LOTS.
6. LINER MATERIAL: 1.02 WHITE POLYETHYLENE FOAM WITH 0.13 NATURAL PTFE-FACED LAMINATED TOGETHER.
7. CAPACITY: 4 OZ, 125.69CC±2.77CC OF

REV	DESCRIPTION OF CHANGE	DRAWN BY	DATE	ECN
B	REVISED	TP	11/02/11	2011101211223320
C	REVISED	TP	02/27/14	2014012714492591
D	REVISED	TP	01/14/16	2015090315204659
-	UPDATE BORDERS & STDS	JD	08/03/22	N/A





FINISH DETAIL
SCALE 2 : 1

ITEM	QTY.	MATERIAL: DESCRIPTION AND SPECIFICATIONS
3	24	50077397 LABEL (NOT SHOWN)
2	24	48-400 CAP W/LINER: WHITE POLYPROPYLENE - SEE NOTE 6
1	24	4 OZ BOTTLE: TYPE III CLEAR SODA-LIME GLASS

THIRD ANGLE PROJECTION

TITLE: 4 OZ BTL, WM, TYPE III CLR, PP/PTFE LNR

THIS CONFIDENTIAL AND PROPRIETARY DRAWING AND THE INFORMATION IN IT ARE THE EXCLUSIVE PROPERTY OF DWK LIFE SCIENCES. IT IS PROVIDED WITH THE UNDERSTANDING AND INTENTION THAT NOTHING CONTAINED HEREIN WILL BE REPRODUCED OR DISCLOSED TO THIRD PARTIES WITHOUT THE PRIOR EXPRESS WRITTEN CONSENT OF DWK LIFE SCIENCES.

ELECTRONIC OR PAPER COPIES OF THIS PRINT OUTSIDE OF DWK LIFE SCIENCES ARE UNCONTROLLED AND ARE FOR REFERENCE ONLY UNLESS OTHERWISE STAMPED.

TOLERANCES
DIMENSIONS IN MILLIMETERS UNLESS SPECIFIED
X ± 2mm
X.X ± 0.5mm
X.XX ± 0.15mm

DURAN WHEATON KIMBLE

ITEM NUMBER:	DRAWING NUMBER:	DRAWN BY:	DRAWN DATE:
W216929	CSD-W216929	BL	05/05/04

SAP NUMBER:	OLD MATERIAL NUMBER:	SHEET SIZE:	SHEET NUMBER:	SHEET SCALE:
W216929	W216929	B	1 OF 1	1:1