Product Specification Sheet

VWR QUALITATIVE FILTER PAPER, GRADE 417

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionés respectent toutes les specifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Catalogue N°: **516-0824**

Application: **Rapid filtrations where fine particle retention is not required.**
Used for coarse particles
Essential oils
Viscous liquids

Material: **100% cellulose with binder**

Technical Specification	Values	Units	Methods
Weight	**110**	**g/m²**	**150536**
Thickness	**creped**	**mm**	
Filtration velocity	**5**	**sec.**	**DIN 53137**
	27	**sec.**	**Herzberg**
Pore size	**35-45**	**µm**	
Wet bursting strength	**15**	**psi**	**interne**
Ash content	**0.09**	**%**	**DIN 53138**

Disclaimer: *VWRwants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there's any doubt, ask us.*