Quantitative Filter Paper, Grade 388
Order Number: FT-*-101-*

Black dot
Fast filtering, wide-pore, lose structure,
ash-free, wet-strengthened

Basis weight [g/m²]:	84
Filtration speed [s]:	10
Wet bursting strength [kPa]:	30
Thickness [mm]:	0.21
Ash content (800°C) [%]:	< 0.01
Chloride content (hot extract) [%]:	< 0.006
pH value (hot extract):	6
Typical retention (*) [µm]:	12-15

(*) Approximate value, not part of the product specification.