Specification

1.50825.0001 LiChrospher® 100 CN (5 µm) LiChroCART® 125-4
HPLC cartridge

	Spec. Values	
Pressure	≤ 30	bar
Theoretical Plates (N/m) (2-Nitrotoluene)	≥ 40000	
Symmetry (2-Nitrotoluene)	0.8 - 1.4	
Capacity factor (2-Nitrotoluene)	0.5 - 1.5	

Peter Schaub

responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature