Technical DataSheet

VWR® ECO VISITOR COAT 30G PP VELCRO CLOSURE ELASTIC WRIST

EU cat. No: 114-9799, 114-9800, 114-9801, 114-9802

CHARACTERISTICS

- Made of 30g non-woven SPP fabric
- White colour
- Size from S to XL
- Nature Rubber Latex Free
- Velcro closure with elastic wrist
- Comply with European Commission Directive 1935/2004/EC, Regulation (EU) No.10/2011, Regulation (EU) No.2016/1416 and Regulation (EU) 2017/752.


PHYSICAL PROPERTIES

Item		Test method	Standard
Weight(g/m2)		ISO3801	30±4
Break Tensile Properties (N/5cm)	MD	ISO13934-1	>40
	CD	ISO13934-1	>25
Elongation Rate(%)	MD	ISO13934-1	>30
	CD	ISO13934-1	>30

Imported by VWR International bv
Address: Researchpark Haasrode 2020, Geldenaaksebaan 464, B-3001 Leuven
E-mail: vwr.be@vwr.com

Ver. 01

PRODUCTS AND DIMENSION SPECIFICATION



Unit: CM

ECN	Colour	Size	A	B	C	D	E	F	G	G'
114-9799	White	S	96	63	45	20	18	48	9	18
114-9800	White	M	100	66	45	22	18	50	9	18
114-9801	White	L	104	69	50	24	18	52	10	20
114-9802	White	XL	108	72	50	26	18	54	10	20

Remark: G: Relax, G': Stretch

PACKAGING AND STORAGE

Package/Quantity	1pcs/bag, 50bags/case
Best Storage Conditions	Temp : 0- 50°C,
	Humidity 45-85%
	without sunshine exposure
Country of Origin	China

QUALITY ASSURANCE CRITERIA

1. Dimension and Function - Critical (AQL1.5, S-2)
- Product dimension out of spec.
- Fabric characteristics out of spec

2. Visual Inspection/ Attributes - Major (AQL 2.5, L-II)
- Foreign Material
- Dirt
- Thread ends loose or too long(over 1cm)
- Holes/Tears
- Foul odor
- Discolor
- Wrong collar tag
- Wrong stitch
- Cloth surface applied wrongly
- Less than 6 stitches per inch
- Dimension tolerance not meet ±1inch for elastic and ±0.5inch for rest

3. Visual Inspection/Attributes - Minor (AQL 4.0, L-II)
- Sewing sundries
- Elastic band end exposed
- Sewing end falling-off
- Jumping Stitch
- Missing stitch
- Unwanted folds
- Elastic tightness non-consistent

Date of Issue: Nov 18, 2025