Identification of the product

Product number	1216652
Product type	Microscope slides with reaction wells
GMDN - Code	57925 General specimen microscopy slide IVD, single-use
GMDN - Definition	A glass or plastic microscopy slide intended to be used for the collection, and preservation and/or transport, of any type of clinical specimen or smear for analysis and/or other investigation. This is a single-use device.
HS-Code	70179000
UDI-DI (smallest sales unit)	04250317302908
UDI-DI (Multipack)	

Manufacturer / Supplier Identification

Manufacturer	Paul Marienfeld GmbH & Co.KG Am Wöllerspfad 4, 97922 Lauda-Königshofen, Germany
Country of origin	Germany

Specification & Compliance

Norm	DIN ISO 8037-1
Dimensions	~76x26 mm
Glass thickness	thickness 1 mm
Material	white glass
Edge processing	edges ground at 90°
Type of marking area	with frosted end
Number of wells x diameter	10 wells of 5 mm diameter
Mask material	white colour mask
Type of packaging inside	in boxes of 50 pcs.
Packing/sales unit	2500 pcs./carton
CE-IVD	CE marked
Brand	MARIENFELD SUPERIOR

Packaging

Smallest sales unit	50 pieces
Multipack	