ReliaPrep™

760501

Characteristics:

Syringe filters with Polyethersulfone (PES) membrane.
Individually packed and sterilized by ethylene oxide.
50 units per pack.

Applications:

Sterile filtration and clarification of aqueous solutions, nutrient media and buffers.

Physical properties	Units	Typical Value
Diameter	mm	28
Filtration Surface	cm^2	6.2
Pore Size	µm	0.2
Housing		MBS- copolymer
Outlet		Male Luer Slip
Inlet		Female Luer Lock

25.06.2019 Rev01

All indicated data to be understood as typical average values.

This product is manufactured in conformance with the Quality Management standard DIN EN ISO 9001.

Contact office:

Ahlstrom-Munksjö Germany GmbH
Niederschlag 1
D-09471 Bärenstein
Germany

T: +49 (0) 37347 / 83-0
F: +49 (0) 37347 / 83-64
Email: filtrak@ahlstrom-munksjo.com

www.ahlstrom-munksjo.com

DISCLAIMER: The information supplied in this document is for guidance only and should not be construed as a warranty. All implied warranties are expressly disclaimed, including without limitation any warranty of merchantability and fitness for use. All users of the material are responsible for assuring that it is suitable for their needs, environment and end use. All data is subject to change as Ahlstrom-Munksjö deems appropriate. Refer to www.ahlstrom-munksjo.com for contact information.

© Ahlstrom-Munksjö 2017