BAKERBOND spe, Alumina

BAKERBOND spe

Disposable Extraction Columns

Material No: 7214-07
Effective from: 2019/10/29

Product Specification

Test	Specification
Color	White
Appearance	Free Flowing Powder
Odour	None (when dry) mineral like (when wet)
pH	7.2 (5% aqueous slurry)
Activity	According to Brockman
Methylene blue exchange	None
Orange GG exchange	None
Bulk Density (g/cc)	0.8 - 0.92
Mean particle size	32 - 63 um
Surface (BET)	200 m2/g
Water Soluble Substances	<= 0.1 %
Purity	99.9 %

For Laboratory, Research or Manufacturing use.
6mL Solid Phase Extraction Columns
1000 mg per Column

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000