Ethanol absolute

AR

Material No: 6777-16
Effective from: 2014/03/20

Product Specification

Test	Specification
Assay (GC)	>= 99.8 %
Water (KF)	<= 0.2 %
Titrable Acid (meq/g)	<= 0.0002
Titrable Base (meq/g)	<= 0.003
Residue after Evaporation	<= 0.001 %
Aldehydes (as CH3CHO)	<= 0.001 %
Higher alcohols (as amyl alcohol)	<= 0.015 %
Methyl Alcohol	<= 0.01 %
Time for reducing KMnO4 solution	>= 10 min
Substances Darkened by H2SO4	passes test
Carbonyl compounds (as CO)	<= 0.003 %

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000