

CATALOG NUMBER	NOMINAL SIZE (mm)	MATERIAL	MEMBRANE MATERIAL
4114-0038	38-430	BLUE HDPE	0.2µ PTFE
4114-0045	45-430	BLUE HDPE	0.2µ PTFE

THREAD FINISH (mm).:	38-430		45-430	
THREAD (T.P.I.)....:	6		6	
	INCH	MM	INCH	MM
A (DIA).....:	1.66	42.2	1.96	49.8
B (HEIGHT).....:	1.17	29.7	1.10	27.9
C (DIA).....:	0.50	12.7	0.50	12.7





SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETORY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	APP'D		REV.
PB	BR 7/9/09	PSS 4114-00XX	0709
CK'D BR			

DESCRIPTION:

VENTED CLOSURES