Product name : **White Filter Paper**	Catalogue no.: **M965FW**
Edition date : 2012-01-12 Approved by : Annette Roy, QA coordinator	Revision : 2013-06-06
This document replaces any previous version	

↺ Product identification

- Designation: **Disposable white filter paper for single funnel**

↺ Production Information

- Producer's name: **Simport *Scientific***
- Address: 2588 Bernard Pilon, Beloeil (Québec), J3G 4S5
- Tel. no / web site: (450) 464-1723 / www.Simport.com

↺ Packaging

- Case: One cardboard box containing 200 units

↺ Traceability.

- Lot no
  - ↪ Composition: 8 digits
  - ↪ Location: ① On exterior case label.

↺ Standard conformity.

- **MIL-STD-105E**: Sampling and inspection procedures
- **FDA**: Filter paper conforms to 21CFR 176.170 and 176.180

Product name : **White Filter Paper**	Catalogue no.: **M965FW**

↺Quality assurance

- No presence of contamination in paper
- Visual attributes

↺ Specifications

- Material: Paper grade which is employed in the foodstuffs sector
- Capillary Rise :120mm|10min
- Differential Pressure at 75 cm/s : 23,0

↺ Storage conditions

- Store in cool dry place.
- Avoid temperature variations and humidity.
- Protect from any possible contamination.

↺ Recommended use

- Recommended for use in the Shandon Cytospin 4 Cytocentrifuge.
  For use with Double CytoSep Cytology Funnel

***OWNERSHIP OF MATERIALS**:*
Materials and information contained in this document are copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.