Chloroform

AR

Stabilized with about 1% Ethanol

Material No: 4444-45
Effective from: 2013/04/12

Product Specification

Test	Specification
Assay	>= 99.5 %
Water (KF)	<= 0.05 %
Titrable Acid (meq/g)	<= 0.0005
Residue after Evaporation	<= 0.001 %
Aldehydes and ketones (as C3H6O)	<= 0.005 %
Substances Darkened by H2SO4	passes test
Free Chlorine (Cl2)	<= 0.0005 %
Stabiliser (Ethanol)	~1 %

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000